

KERN 326-12

Single weight, compact shape, polished stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	1 x 1 x 1 mm
Gross weight:	2 kg
Net weight:	2 kg
Shipping method:	Parcel service

Construction

Material:	polished stainless steel
Shape:	compact shape

Measuring system

Nominal value:	2 kg
OIML class:	F1
Tolerance:	10 mg

Accessories

Plastic box:	347-120-400
Wooden box:	317-120-100

Service

DAkkS Certificate:	962-442
Verification:	952-442