

Clifton Unstirred Water Baths

NE1D Series Unstirred Digital Water Baths

Clifton Range®



NE1D-14



Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. The NE1D series, Clifton Digital Baths are designed for immersion of flasks/bottles/racks in a highly stable temperature environment throughout the water bath.

NE1D-4

NE1D-8

NE1D-14

NE1D-22

	NE1D-4	NE1D-8	NE1D-14	NE1D-22
GENERAL				
Temp. Range	Ambient + 5°C - 99°C			
Temp. Control	Digital, LED display showing actual or set temperature, resolution 1°C			
Sensitivity/Uniformity	±0.2°C / ±0.1°C			
Heater Power (220-240V)	400W	800W	1000W	1500W
Voltage	120V or 230V models available 50-60Hz			
Drain (12mm bore)				Drain Outlet
DIMENSIONS				
Tank Capacity	4 Litre	8 Litre	14 Litre	22 Litre
External (mm)	332w x 185d x 290h	332w x 270d x 290h	361w x 332d x 290h	538w x 332d x 290h
Internal (mm)	300w x 150d x 150h	300w x 240d x 150h	325w x 300d x 150h	500w x 300d x 150h
Water Level Max. (mm)	130			
Water Level Min. (mm)	55			

NE1D Series Unstirred Digital Water Baths

Water Baths Accessories - Lids



Stainless Steel Gabled lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	SL1-4	SL1-8	SL1-14	SL1-22



Stainless Steel Flat One Piece lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	LD-4	LD-8	LD-14	LD-22



Stainless Steel Hinged Gabled lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	N/A	N/A	SL1-14H	SL1-22H



Clear Hinged Gable Lid

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	GL1-4	GL1-8	GL1-14	GL1-22



Clear Lift Off Gable Lid

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	GL1-4AF	GL1-8AF	GL1-14AF	GL1-22AF



Stainless Steel Lid with Concentric rings

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	SLR1-4	SLR1-8 SLR2-8	SLR1-14 SLR2-14 SLR3-14	SLR1-22 SLR3-22

SLR1-4: 2 x 105mm holes

SLR1-8: 4 x 83mm holes

SLR2-8: 2 x 105mm holes

SLR1-14: 4 x 105mm holes

SLR2-14: 6 x 83mm holes

SLR3-14: 4 x 114.5mm holes

SLR1-22: 5 x 105mm holes

SLR3-22: 6 x 114.5mm holes



Polypropylene Spheres (Pack of 200)

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
No. of packs	1 Pack	2 Pack	2 Packs	4 Packs
Accessory Code	BP0368			