

# Cole-Parmer®

## Cole-Parmer® Unheated Baths - WBU-200-8, WBU-200-12, WBU-200-18, WBU-200-26, WBU-200-48

- Stainless steel construction
- Seam-free and corrosion resistant stainless steel inners for easy cleaning
- Rugged splash-proof case
- Integrated carrying handle
- Maximum working temperature of 200°C
- All models come with a 3 year warranty as standard



Unheated Bath WBU-200-8

## Cole-Parmer® Unheated Baths - WBU-200-8, WBU-200-12, WBU-200-18, WBU-200-26, WBU-200-48

These Cole-Parmer® are designed to be used with a clip-on Tempette or Tempunit® thermoregulator, these water baths incorporate carrying handles for added safety. All water baths have stoved enamelled steel outer cases and are supplied with bridge mounting plate to hold the thermoregulator.



Unheated Bath WBU-200-8



Unheated Bath WBU-200-12



Unheated Bath WBU-200-18



Unheated Bath WBU-200-26



Unheated Bath WBU-200-48

## Technical Specification

Specification	WBU-200-8	WBU-200-12	WBU-200-18	WBU-200-26	WBU-200-48
Capacity litres	8	12	18	26	48
Dimensions (mm) Length	265	354	530	530	594
Width (mm)	325	325	325	325	365
Height (mm)	172	172	172	222	298
Internal Dimensions Length (mm)	240	329	505	505	559
Width (mm)	300	300	300	300	330
Height (mm)	150	150	150	200	274
Top of bath to liquid level max depth (mm)	65	65	65	65	65
Working length to thermoregulator (mm)	115	205	380	380	430
Working depth - max/min (mm)	130/100	130/100	130/100	180/150	255/224
Working capacity - max/min (litres)	8.0/6.0	11.6/8.4	18.0/13.2	26.0/20.5	48.5/42.5
Shipping Weight	5.5kg	6.1kg	7.5kg	9.5kg	14.6kg

## Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® Stainless steel bath, 8 litre capacity	16104-53	WBU-200	WBU-200-8	FBATH08
Cole-Parmer® Stainless steel bath, 12 litre capacity	16101-88	WBU-200	WBU-200-12	FBATH12
Cole-Parmer® Stainless steel bath, 18 litre capacity	16112-01	WBU-200	WBU-200-18	FBATH18
Cole-Parmer® Stainless steel bath, 26 litre capacity	16101-89	WBU-200	WBU-200-26	FBATH26
Cole-Parmer® Stainless steel bath, 48 litre capacity	16104-45	WBU-200	WBU-200-48	FBATH48

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials  
an antylia scientific company

Find out more at antylia.com