

IKA

designed for scientists



DB 1.3

/// Data Sheet

Used for VWR microtubes (1.5 ml)

Pore size: Ø 10.9 mm

No. of holes: 20

Hole depth: 39.1 mm

Dimensions (W x D x H): 95 x 76 x 51 mm

www.ika.com

Subject to technical changes





designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.858
Permissible ambient temperature [°C]	5 - 40