

IKA

designed for scientists



DB 2.2

/// Data Sheet

Used for Conical tubes (50 ml)

Pore size: Ø 29.0 mm

No. of holes: 5

Hole depth: 47.6 mm

Dimensions (W x D x H): 95 x 76 x 51 mm



designed for scientists

Technical Data

| | |
|--------------------------------------|--------------|
| Dimensions (W x H x D) [mm] | 95 x 51 x 76 |
| Weight [kg] | 0.6298 |
| Permissible ambient temperature [°C] | 5 - 40 |