

# IKA

designed for scientists



## DB 4.5

/// Data Sheet

Used for round bottom tubes (15/16 mm)  
Pore size: Ø 17.5 mm  
No. of holes: 12  
Hole depth: 48.4 mm  
Dimensions (W x D x H): 95 x 76 x 51 mm

[www.ika.com](http://www.ika.com)

Subject to technical changes





designed for scientists

### Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.651
Permissible ambient temperature [°C]	5 - 40