Reducing Patient Identification Errors

while enhancing lean workflow

Thermo Scientific PrintMate AS Cassette Printer

Eliminate waste and streamline your labeling process by printing cassettes on-demand with the Thermo Scientific[™] PrintMate[™] AS Cassette Printer at the grossing station.

- Prints each cassette in 10 seconds or less
- Small footprint allows the printer to fit easily next to the grossing station
- Print only the cassettes you need, right at the point of use
- Print barcodes on cassettes to automate downstream processes
- Eliminate manual handwriting processes





Are specimen identification errors a concern in your laboratory?

The small footprint of the Thermo Scientific™ SlideMate™ AS Slide Printer brings the convenience of on-demand slide printing directly to the microtome. The SlideMate AS Slide Printer drastically reduces the chance for patient identification errors by allowing the laboratory to print accurate information at the moment the slide is required.





Order Information		
Description	Order Number	
PrintMate AS 150 Cassette Printer	A83800001	
PrintMate AS 450 Cassette Printer	A83800002	
PrintMate AS 900 Cassette Printer	A83800003	
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), LabWriter software CD, operator guide		
PrintMate AS 150 Cassette Printer	A83800004	
PrintMate AS 450 Cassette Printer	A83800005	
PrintMate AS 900 Cassette Printer	A83800006	
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), Thermo Scientific™ LabWriter™ software CD, operator guide, all-in-one PC, digital imager		

PrintMate AS Accessories	
Description	Order Number
Collection system (includes cable, 4 trays for cassettes, output guides)	A79510014
Collection trays (box of 4)	A79510022
Hot foil tape, box of 5 rolls	A79510013
All-in-one touch screen PC	B81300033
Digital imager (linear/2D barcode reader)	B35310050
Cassette hopper (ea)	A79510024

Specifications Specification Specif		
Dimensions: (D x W x H):		
PrintMate AS 900 w/o cassette t	ubes 15 x 17 x 11.4 in. (38 x 43.5 x 29 cm)	
PrintMate AS 900 w/ cassette tu	bes 15 x 17 x 32 in. (38 x 43.5 x 81 cm)	
Weight:		
PrintMate AS 900 w/o cassette t	ubes 31 lb (14 kg)	
PrintMate AS 900 w/ cassette tu	bes 37.5 lb (17 kg)	
Power Requirements:	100–240 VAC, 50/60 Hz 150 VA (512 BTU/h)	
Print media type:	Thermal transfer ribbon	
Print media capacity:	Approximately 20,000 cassettes per rol	
Print Speed:	10 seconds per cassette	
Connectivity:	USE	
Print media color:	Black	
Cassette Capacity:		
PrintMate AS 150	150 cassettes	
PrintMate AS 450	450 cassettes	
PrintMate AS 900	900 cassettes	

To learn more about how on-demand cassette printing can impact in your laboratory, call: (800) 522-7270

www.thermoscientific.com/pathology

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Anatomical Pathology

4481 Campus Drive Kalamazoo, MI 49008 UNITED STATES

Tudor Road, Manor Park Runcorn, WA7 1TA UNITED KINGDOM +1 (800) 522-7270

+44 (0) 800 018 9396 +44 (0) 1928 534 000

