















Overview: Supported inkjet printers (1)

Discontinued, supported printers: MI 9020, MI 9232	MI 9028	MI 9410	MI 9450
	 www.markem-imaie.com/	 www.markem-imaie.com/	 www.markem-imaie.com/
Performance	<ul style="list-style-type: none"> ▪ 4 Lines ▪ Alphanumeric ▪ Graphics/Logo printing ▪ Side wall printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing
Other features	<ul style="list-style-type: none"> ▪ Bench top ▪ Key pad ▪ USB port ▪ Ethernet 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB and SD port ▪ Ethernet 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB and SD port ▪ Ethernet
Dimensions (HxWxD)	595 x 307 x ??? mm	555 x 360 x 357 mm	555 x 360 x 357 mm
Weight (kg)	18.5	25	25
List Price (EUR)	<10'000	>10'000	>10'000
Maintenance	<ul style="list-style-type: none"> ▪ Maintenance: every 24 month ▪ Head cleaning recommended every second week ▪ Head can be left attached on MEDIAJET 	<ul style="list-style-type: none"> ▪ Maintenance: every 18 month ▪ Head cleaning recommended every second week ▪ Head can be left attached on MEDIAJET 	<ul style="list-style-type: none"> ▪ Maintenance: every 18 month ▪ Head cleaning recommended every second week ▪ Head can be left attached on MEDIAJET
INTEGRA comments	Industrial interface must be ordered.	Recommended MI printer , lower consumable and maintenance costs. Industrial interface must be ordered.	Industrial interface must be ordered. Higher line speed and storage capacity.

Overview: Supported inkjet printers (2)

Discontinued, supported printers: Linx 4900, Linx 5900, Linx CJ400 Domino A320i, Domino A420i	Linx 8900 (89xx: -10, -20, -40)	Domino Ax150i	Domino Ax350i
	 www.linxglobal.com	 www.domino-printing.com	 www.domino-printing.com
Performance	<ul style="list-style-type: none"> ▪ 8 Lines ▪ Alphanumeric, barcodes, 2D codes ▪ Graphics/Logo 25 px ▪ Side wall / bottom printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing
Other features	<ul style="list-style-type: none"> ▪ Benchtop ▪ Touch screen ▪ USB port ▪ Ethernet 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB port ▪ Ethernet, RS232, GPIO 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB port ▪ Ethernet, RS232, GPIO
Dimensions (HxWxD)	503 x 440 x 355 mm	373 x 400 x 393 mm	411 x 430 x 331 mm
Weight (kg)	24	16	22
List Price (EUR)	>10'000	8'300	>10'000
Maintenance	<ul style="list-style-type: none"> ▪ Maintenance: every 24 months ▪ Head cleaning every 3 months 	<ul style="list-style-type: none"> ▪ Low-maintenance ▪ Head can be left attached on MEDIAJET 	<ul style="list-style-type: none"> ▪ Low-maintenance ▪ Head can be left attached on MEDIAJET
INTEGRA comments	This reliable Linx printer replaces all previous models.	Recommended Domino printer Recommended, if bar codes or 2D codes needed.	Recommended, if bar codes or 2D codes needed.

Ordering information MEDIAJET inkjet printer connections kits

Discontinued, supported printers: MI 9020, MI 9232, Linx 4900, Linx 5900, Linx CJ400 Domino A320i, Domino A420i	Markem-Imaje (MI) 9020/9028	Markem-Imaje (MI) 9232/9410/9450	Linx 4900/5900/8900/ CJ400 (side wall printing)	Linx 4900/5900/8900 (bottom printing)	Domino A320i/A420i/ Ax150i/Ax350i
					
	Part number				
Printer connection kit: with printer head fixing device and interface cable 	103080	103091	113845	113840	113060
Cart 	-	-	-	113841 (Linx 4900/5900 only)	-
Interface cable (spare part) 	103081	103081	113853	113853	103082
INTEGRA comments	Compatible with Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes.	Compatible with Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes. Connection kit does not fit to Linx CJ400.	Compatible with Ø 90 mm or Ø 60 mm Petri dishes.