

BradyJet J2000 Colour Label Printer

Identify samples by colour to increase productivity



www.bradyeurope.com/J2000











The more samples laboratories need to process, the more time consuming it becomes to retrieve one specific sample. Brady Corporation proposes sample colour coding to solve this and offers the BradyJet J2000 Colour Label Printer to do it, to stay productive and to avoid spending too much time on sample hunting.

Sort samples by colour

The BradyJet J2000 Colour Label Printer can print any colour on reliable laboratory sample labels. By setting up a colour code, laboratory specialists will be able to recognise sample categories at a glance. With over 16 million colour variations, users may switch colour every day, or use a different colour per application, storage or urgency type depending on the needs of their laboratory. Next to colour, the BradyJet J2000 also offers 4800 dpi photo quality print to add highly legible data or barcodes a decent scanner could never misread.

Print reliable labels

The BradyJet J2000 prints on Brady's reliable laboratory sample labels that come in matching sizes and shapes for most sample containers including tubes, vials, cryo containers, dishes, plates and general glassware. These labels are designed to stay attached and to remain legible when printed with a Brady printer, even if they are stored in liquid nitrogen and in freezers, or when used in hot water baths or autoclaves. Brady's laboratory sample labels can also resist common laboratory chemicals like acetone, ethanol, toluene, xylene and IPA.

Easy to use

With 26.4 by 38.9 cm dimensions, the BradyJet J2000 has a small footprint and is very practical and easy to use. Spring loaded Brady label rolls can be switched in less than 20 seconds and with automatic label set up, the printer is ready to print in no time. And with up to 63.5 mm per second print speed, even larger colour coded sample labels can quickly be printed on demand.



For more information or a live demo, visit

www.bradyeurope.com/J2000

Printer specifications

PRINTER				
Description	BradyJet J2000 Colour Label Printer			
Includes	This kit contains a J2000 printer, a power cord and AC adapter, a USB cable, a Quick Start Guide, an ink cartridge, a material roll for printhead alignment and a Brady Workstation Suite			
Print Speed (mm/sec)	63.5			
Print Technology	Inkjet			
Print Resolution (dpi)	Up to 4800			
Stand-alone Printing (Detached from PC)	Peripheral			
Colour Capability	Full Colour			
Weight (kg)	3.20			
Size W x H x D (mm)	264 x 180 x 389			
Cutter	Auto Cutter			
Warranty	1 year			
MEDIA				
Continuous vs Die-cut Labels	Continuous Tape & Die Cut Labels			
Custom Labels	Yes			
Materials Supported	B-2595, B-2569, B-2551, B-7425J			
Min. Label Length (mm)	6.35			
Min. Label Width (mm)	19			
Max. Label Width (mm)	101.60			
Max. Liner Width (mm)	108.00			
Max. Print Width (mm)	101.60			
Recommended Usage per Day	Around 3000 labels per day, depending on specific label size, text printed and other factors			
SOFTWARE & CONNEC	TIVITY			
Connectivity Options	USB			
Software Compatibility	Brady Workstation, Windows based driver for 3rd party software use			
COMPLIANCE				
Approval/Compliance	USA: UL, FCC Class A; Canada: cUL, IC/ISED; European Union : CE, RoHS 2011/65/EU, 2014/35/EU, LVD, 2014/30/EU, RED 2014/53/EU; Australia/			

Printers

Order Ref.	Description
J2000-EU	BradyJet J2000 Colour Label Printer - EU
J2000-EU-LABS	BradyJet J2000 Colour Label Printer - EU with Brady Workstation LAB Suite
J2000-UK	BradyJet J2000 Colour Label Printer - UK
J2000-UK-LABS	BradyJet J2000 Colour Label Printer - UK with Brady Workstation LAB Suite
J2000	BradyJet J2000 Colour Label Printer - US
J2000-US-LABS	BradyJet J2000 Colour Label Printer - US with Brady Workstation LAB Suite
	J2000-EU J2000-EU-LABS J2000-UK J2000-UK-LABS J2000



The Basic Suite is included with every printer, all you need to do is download it. Workstation.BradyID.com/Free



Easily identify any sample with the most wanted laboratory sample label design apps in our Brady Workstation Lab suite.

Inks

Art. No.	Order Reference	Description
150154	J20-CMY	CMY pigment based ink cartridge for J2000 Printer

Accessories

Art. No.	Order Reference	Description
150160	J20-CUTTER	Replacement Cutter Blade for J2000
150163	J20-PAD	Replacement Ink Pads for J2000





For more information on supported LAB materials, visit www.bradyeurope.com/J2000

Materials

Diagram	Order Reference	Art. No.	Colour	Label Width A (mm)	Label Height B (mm)	Diameter C (mm)	Total Qty/ Pack	
Inkjet Receptive Polypropylene								
A	J20-257-7425J	150007	White	20.96	6.35	_	1500	
A	J20-151-7425J	150011	White	50.80	6.35	_	2000	
A	J20-152-7425J	150005	White	25.40	9.53	-	2000	
•	J20-249-7425J	150023	White	25.40	9.53	9.53	1500	
A	J20-59-7425J	150006	White	25.40	12.70	-	2000	
-	J20-251-7425J	150033	White	25.40	12.70	11.18	1500	
A	J20-37-7425J	150020	White	38.10	12.70	-	1500	
_	J20-53-7425J	150024	White	50.80	12.70	-	2000	
A	J20-256-7425J	150016	White	33.02	15.24	-	1500	
A	J20-88-7425J	150015	White	25.40	19.05	-	2000	
A	J20-136-7425J	150019	White	38.10	19.05	_	1500	
A	J20-179-7425J	150025	White	25.40	25.40	-	1500	
A	J20-235-7425J	150022	White	38.10	25.40	-	1250	
	J20-137-7425J	150018	White	50.80	25.40	_	1500	
A	J20-258-7425J	150028	White	38.10	44.45	-	900	



EUR-M-001-EN 10/06/2019 © 2019 Brady Worldwide, Inc. All Rights Reserved.