

# PATHOSHIELD



**SAMPLE  
TRANSPORT  
PACKAGING**



The PathoShield range are complete shipping solutions designed to carry a wide range of primaries, including blood vials, swabs and sample containers. Each system is pre-printed with the required markings and includes the relevant components to ship samples in compliance with transport regulations.

## THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

## FEATURES

- Complete Category A and Category B solutions include flexible 95kPa tested compliant secondary container, absorbent, labels, and UN markings
- Range of sizes suitable for shipping through postal and courier networks
- Carry a wide range of primaries including blood vials, swabs and sample containers
- Can be supplied as an insulated pack
- Supplied flat packed for efficient storage and shipping prior to use

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



## PATHOSHIELD RANGE

Product Code	Product Name	External Dimensions	Pouch Internal Dimensions (W x L)	Number of Pouches	Number of Tubes	Mass	Volumetric Weight L/6
--------------	--------------	---------------------	-----------------------------------	-------------------	-----------------	------	-----------------------

### Category A Solutions

PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	3	0.23 kg 0.5 lbs	0.9 kgvol <sup>6</sup> 2.1 lbsvol <sup>6</sup>
-------	--------------------------	---	------------------------------	---	---	--------------------	---

### Category B Solutions

PA003	PathoShield 3 w A4 Pouch	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	6	0.16 kg 0.35 lbs	0.5 kgvol <sup>6</sup> 1.0 lbsvol <sup>6</sup>
PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	6	0.23 kg 0.5 lbs	0.9 kgvol <sup>6</sup> 2.1 lbsvol <sup>6</sup>
PA001	PathoShield 5 w A5 Pouch	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.09 kg 0.20 lbs	0.4 kgvol <sup>6</sup> 0.8 lbsvol <sup>6</sup>
PA006	PathoShield 6 w A5 Pouch	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.07 kg 0.15 lbs	0.2 kgvol <sup>6</sup> 0.4 lbsvol <sup>6</sup>
PA007	PathoShield 7 w A5 Pouch	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	155 x 180mm 6.1 x 7.0 in	1	4	0.10 kg 0.21 lbs	0.3 kgvol <sup>6</sup> 0.6 lbsvol <sup>6</sup>
PA008	PathoShield 8 w A5 Pouch	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	155 x 180mm 6.1 x 7.0 in	1	4	0.04 kg 0.09 lbs	0.1 kgvol <sup>6</sup> 0.2 lbsvol <sup>6</sup>
PA023	PathoShield 3 w A4 Bag	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.13 kg 0.28 lbs	0.5 kgvol <sup>6</sup> 1.0 lbsvol <sup>6</sup>
PA024	PathoShield 4 w A4 Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.19 kg 0.43 lbs	0.9 kgvol <sup>6</sup> 2.1 lbsvol <sup>6</sup>
PA025	PathoShield 5 w A5 Bag	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.08 kg 0.18 lbs	0.4 kgvol <sup>6</sup> 0.8 lbsvol <sup>6</sup>
PA026	PathoShield 6 w A5 Bag	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.06 kg 0.13 lbs	0.2 kgvol <sup>6</sup> 0.4 lbsvol <sup>6</sup>
PA027	PathoShield 7 w A5 Bag	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.09 kg 0.19 lbs	0.3 kgvol <sup>6</sup> 0.6 lbsvol <sup>6</sup>
PA028	PathoShield 8 w A5 Bag	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.04 kg 0.09 lbs	0.1 kgvol <sup>6</sup> 0.2 lbsvol <sup>6</sup>

### Outer

PA014	PathoShield 4 Outer Box	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	-	-	-	0.16 kg 0.4 lbs	0.9 kgvol <sup>6</sup> 2.1 lbsvol <sup>6</sup>
PA015	PathoShield 5 Outer Box	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	-	-	-	0.07 kg 0.15 lbs	0.4 kgvol <sup>6</sup> 0.8 lbsvol <sup>6</sup>

### Mailing Bag

MBUN3373C5*	Mailing bag UN3373 C5 1000 Units	240 x 170 mm 9.4 x 6.6 in	-	-	-	-	-
-------------	-------------------------------------	------------------------------	---	---	---	---	---

\*UN3373 Mailing Bag is designed to be used in conjunction with a fully compliant UN3373 Category B packaging solution

**UNITED KINGDOM:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom

t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICA:** DGP Intelsius LLC  
1399 Perry Rd, Suite 190,  
Plainfield, IN 46168,  
United States of America

t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**MALAYSIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia

t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com

**GERMANY:** DGP Intelsius GmbH  
Eulerweg 11  
64291 Darmstadt  
Germany

t: +49 6151 4936650  
f: +49 (0) 6151 4936651  
info@intelsius.de

**IRELAND:** DGP Intelsius Ireland Ltd.  
Unit 8, Coolcaslagh Industrial Estate,  
Coolcaslagh, Killarney  
Co. Kerry, V93 K006, Ireland

t: +353 (0)64 668 61 71  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
Sun Srishti, C 10, Saki Vihar Road,  
Chandivali Powai, Mumbai - 400072  
Maharashtra - India

+91 22 49746536  
compliance@intelsius.com