



Model: CMRTSG125

## 145L Glass Door RTS Cabinet

### Overview

- +15°C to +25°C temperature range
- Core temperature display, represents your vaccine vial through **iCOOL** INTELLIGENT COOLING temperature probe
- External digital controller and temperature display
- Battery back up
- Annual service & Calibration available
- Fitted lock with 2 keys



**1 YEAR PARTS & LABOUR WARRANTY**

**ANNUAL SERVICE & CALIBRATION AVAILABLE**

### Key Features




#### Controller

Controller system. Read, record, reset with the touch of one button. Adjustable temperature setting.



#### Battery back up

Battery back up should there ever be a power outage the battery backup will continue to monitor the internal temperature of the fridge.



#### Porthole

Factory fitted porthole enables quick and easy access for any secondary temperature monitoring probes.



#### LED lighting

Low voltage LED lighting improves running cost as no heat is produced. Long lasting.



#### Flexible storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.


## Details

Model: CMRTSG125

|                                                                                                       |                                          |
|-------------------------------------------------------------------------------------------------------|------------------------------------------|
| Gross Capacity (litres)                                                                               | 148                                      |
| Net Capacity (litres)                                                                                 | 145                                      |
| Unit External Width (mm)                                                                              | 595                                      |
| Unit External Depth (mm)                                                                              | 595                                      |
| Unit External Height (mm)                                                                             | 815                                      |
| Product Weight (kg)                                                                                   | 44                                       |
| Packed Width (mm)                                                                                     | 650                                      |
| Packed Depth (mm)                                                                                     | 650                                      |
| Packed Height (mm)                                                                                    | 865                                      |
| Packed Weight (kg)                                                                                    | 47                                       |
| Temperature Range                                                                                     | +15°C to +25°C                           |
| External Finish, Internal Finish                                                                      | Painted Steel, Painted Steel             |
| Battery Back Up                                                                                       | Yes                                      |
| Door Finish, Door Type, Door Style                                                                    | Painted Steel, Glass, Hinged, Right Hand |
| Open Door Sensor                                                                                      | Yes                                      |
| Lock                                                                                                  | Yes                                      |
| Adjustable Feet (number)                                                                              | Yes (2)                                  |
| Porthole                                                                                              | Yes                                      |
| Porthole Location                                                                                     | Rear, Middle                             |
| Temperature Display on Panel, Temperature Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms | Yes                                      |

|                                   |                                     |
|-----------------------------------|-------------------------------------|
| Number of Shelves                 | 4 (2 full size, 1 small base shelf) |
| Number of Shelf Positions         | 6                                   |
| Shelf Material                    | Plastic Coated Wire                 |
| Interior Light                    | Yes                                 |
| Light Type                        | LED                                 |
| Light Activation                  | Switch                              |
| Fan                               | Yes                                 |
| Spark Free                        | No                                  |
| Controller Type                   | Digital, External                   |
| Number of Probes                  | 2                                   |
| Probe Type                        | Air & Load                          |
| Ambient Temperature Range         | +10°C to +30°C                      |
| Refrigerant Type                  | R600a                               |
| CFC Free, HCFC Free, Ammonia Free | Yes                                 |
| Voltage Supply                    | 220-240V/ 50Hz                      |
| Certification                     | CE/UK                               |
| Noise Level (dBA)                 | 50                                  |
| Input Power (W)                   | 76                                  |
| Cable Length (cm)                 | 20                                  |
| Energy Consumption (kWh/annum)    | 0.9                                 |

## Related Products

**CMRTSS125**145L Solid Door  
RTS Cabinet**CMDL01**Advanced Data Logger  
with LCD Display**COMBI01**Keyless  
Combination Lock