

Underbench Sparkfree Refrigerator Model: BIO130FRSS

The *BioCOLD standard* series of refrigeration has been built to a very high specification, to meet the requirement for general laboratory storage.

This Underbench refrigerator features a spark free interior as standard for the safe storage of volatile substances or samples. Options on this model include a door mounted temperature display/alarm module.



Standard Features

- Sparkfree interior
- Reversible doors
- Replaceable door gasket
- CFC and HCFC free refrigeration system

Options

- Lock
- Door mounted temperature display/alarm module

Technical Specification

Dimensions Temperature External (H x W x D), (mm) 845 x 550 x 580 Set Point (adjustable) (°C) +4 Weight (kg) Operating range (°C) +2 to +8 Capacity **Electrical Supply** Refrigerator Volume (L)(Net) 130 220-240 Voltage (V) Shelves (Safety glass) Frequency (Hz) 50 Fuse (A) 13 **Defrost Type** Manual



Telephone/Fax: 01506 873795 E-mail: sales@biocold.co.uk www.biocold.co.uk