

Upright Refrigerator Model: BIO360FRS

The <u>BioCOLD</u> <u>STANDARD</u> series of refrigeration has been built to a very high specification, to meet the requirement for general laboratory storage.

The Upright Refrigerator features 5 glass safety shelves as standard.



Standard Features

- Glass shelves
- Reversible doors
- Replaceable door gasket
- CFC and HCFC free refrigeration system

Options

- Sparkfree interior
- Door mounted temperature display/alarm module

Technical Specification

Temperature Dimensions External (H x W x D), (mm) 1750 x 600 x 655 Set Point (adjustable) (°C) Weight (kg) Operating range (°C) +2 to +8 Capacity Refrigerator Volume (L) 350 **Electrical Supply** 220-240 Voltage (V) Shelves (Glass) 5 Frequency (Hz) 50 Fuse (A) **Defrost Type** Manual

BioCOLD

Telephone/Fax: 01506 873795 E-mail: sales@biocold.co.uk www.biocold.co.uk