





This **Labcold** Sparkfree fridge freezer is engineered to store flammable materials safely with all potential sources of ignition, such as the thermostat removed from inside the compartments of both the fridge and the freezer.

Warning: Do not store open containers of volatile substances in this fridge freezer

Did you know that 25% of chemical accidents are caused by improper storage?



Source: US Chemical Safety and Laboratory Investigation Board

compartment as standard, this fridge freezer can also be factory fitted with a Labcold digital thermometer which has a high and low alarm and records minimum and maximum temperatures.

As you would expect, your sparkfree fridge freezer is backed by a comprehensive UK parts and labour warranty and as Labcold are an accredited UKAS temperature calibration laboratory #8898 we can take care of all your calibration needs as part of a general maintenance contract. Please ask for a quote when you place your order.

RLFF13248LK Sparkfree Fridge Freezer

Part Number	RLFF13248LK
Capacity refrigerator	271 litres
Capacity freezer	113 litres
Style	Freestanding
Exterior Dimensions (HxWxD)	2000 x 625 x 680mm
Weight	67.5kg
Energy consumption	0.88kWh/24hr
Operating temperature	Fridge +0°C to +10°C Freezer -18°C to -23°C
Interior construction	Moulded Plastic
Exterior construction	White Enamel
Door style	Solid
Door hinging	Right
Door alarm	\checkmark Refrigerator only
Shelves	5
Drawers	3
Feet	4
Door lock with 2 keys Refrigerator Freezer	\checkmark
Digital thermometer with temperature display, high and low alarms and min/max recording	Optional
High Temperature Alarm	 ✓ Visual only and only for the refrigerator compartment
Defrost	Fridge - auto Freezer - manual
Low GWP refrigerant	\checkmark
1 year UK parts and labour warranty	\checkmark





UKAS

01256 705 570

www.labcold.com