



# ASAP INTERNATIONAL SDN. BHD. (1157980-X)

GST No : 000696299520

No. 1, Jalan Sitar 33/6, Seksyen 33,  
40400 Shah Alam, Selangor, Malaysia.

Tel: +6 03 5191 0166 Fax: +6 03 5191 0702 Website: whyasap.com

## PRODUCT SPECIFICATION

Product	: ASAP X-Tra Thick X-Tended Black Powder Free Nitrile Examination Gloves
Product Description	: Non-Sterile, Powder Free Nitrile Exam Gloves, Palm Textured, Black Colour, Beaded Cuff, Ambidextrous, For Single Use
Available Sizes	: Extra Small, Small, Medium, Large, Extra Large, Double Extra Large
Packing Mode	: 100 pcs per box, 10 boxes per carton
Quality Assurance	: Meets EN 455 Standards

Test	Particular	Specification	
Dimension (mm)	Palm Width	Extra Small	75 ± 5
		Small	85 ± 5
		Medium	95 ± 5
		Large	105 ± 5
		Extra Large	115 ± 5
		Double Extra Large	125 ± 5
	Length		Min. 290
Thickness (Single wall)	Cuff	min. 0.09	
	Palm	min. 0.12	
	Finger	min. 0.13	
Physical Properties	Force at Break (N)	Before Aging	Min. 6
		After Aging	Min. 6
Physical Appearance	Rolled beaded cuff Black Colour Palm Textured Surface Ambidextrous		

*\*\* Remark: The above data represent the norm of the product which may be changed without notice.*

### Handling and Storage

Store in a dry place at room temperature. Exposed product should be shielded from direct sunlight, intense artificial light, x-ray machines and other sources of ozone.

### Caution

This product contains nitrile rubber which may cause an allergic reaction in some individuals. Do not expose these gloves to any person with a known or suspected sensitivity to nitrile rubber. If allergic reaction occurs, discontinue use immediately and consult a physician. This product is not made with natural rubber latex.

