


GN830

GN830

Indigo nitrile powder free examination glove with a long cuff



Sectors:

- Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet and Dental

Features

Manufactured from: Nitrile butadiene

Accelerator Free: Removing the risk of contact dermatitis caused by these chemicals

Latex free: Reducing the risk of hypersensitivity and allergic reactions

Suitable for use with cytotoxic/chemotherapy drugs

Food Contact: Tested in accordance with European standard EN1186 making it suitable for food contact

Powder free: Reduces the risk of dust contamination

Textured finish: For improved grip

Long Cuff: 29cm for extra protection

Code	GN830
Colour	Indigo
Sizes	S-XL



Inners dimensions (mm)
260 x 126 x 70mm

Outers dimensions (mm)
385 x 260 x 270mm

Size	Inner Barcode	Outer Barcode
S/7	5024951017073	502495107110
M/8	5024951017080	5024951017217
L/9	5024951017097	5024951017134
XL/10	5024951017103	5024951017141

Standards

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with the Medical Devices Regulation (EU) 2017/745 Class 1

Tested to European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Product Food Approved to EC1935/2004

EU Type Examination carried out by SATRA

Tested to European Standards EN420:2003+A1:2009 & EN ISO 374-1:2016 & EN ISO 374-5:2016

Packing

100 gloves per dispenser
10 dispensers per case
40 cases per pallet



MDR CLASS I
PPE CAT III
2777



K



VIRUS



For further technical information, please contact Customer Services on 033 3320 8550 or the Technical team on 033 3320 8555

