


1. Part Number: 1030DT all sizes

2. Product Information:



| Part No: | Nominal Volume | Weight |
|------------|----------------|--------|
| 1030DT-20 | 20ml | 486g |
| 1030DT-25 | 25ml | 456g |
| 1030DT-50 | 50ml | 327g |
| 1030DT-100 | 100ml | 550g |

Description

Tips for use with 32mm diameter PowderThief Bodies

Materials of Construction:

316L stainless steel

Finish:

<1 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

3. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jan-21