


| | | | |
|---|------------------------------|----------------------|----------|
|  | Product Specification | Specification Number | 1320A-01 |
| | | Revision | 01 |

Part No. 1320A-01
Description Sack Master – with Plain End
Material of Construction: 316 stainless steel



Method of Construction: Single piece construction
Surface Finish: Better than 0.5 microns Ra
Overall Length: 400mm
Outer Diameter: 25mm
Inner Diameter: 22mm
Nominal Weight: 335g
BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are of vegetable origin

| | | |
|----------------------------|--|--|
| Specification Approved By: | Signature:  | Name: Jonathan Fenner Position: Director Date: 24-Jan-13 |
| Page 1 of 1 | Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com | |