

Product Specification

Specification Number

1320C-01

Revision

ion 01

Part No. 1320C-01

Description Sack Master – with Bottle Adapter

Material of Construction: 316 stainless steel





Method of Construction:

Bottle adapter screws on to the sack sampler tube (can be easily

unscrewed for cleaning)

Surface Finish: Better than 0.5 microns Ra

Overall Length: 420mm

Outer Tube Diameter: 25mm

Inner Diameter: 22mm

Nominal Weight: 460g

BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are

of vegetable origin

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director
Date: 24-Jan-13

Sampling Systems Ltd.