

1. Part Number: 8060H-1005S

2. Photograph:



3. Product Information:

Description:	SteriWare PowderThief - 1000mm – 5ml Tip - Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	987mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	5ml
Width of Handle:	103mm
Nominal Weight of Sampler:	165g
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25 kGy)
Number of Samplers per Box:	20
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21