

**1. Part Number:** 8060H-1520

**2. Photograph:**



**3. Product Information:**

Description: SteriWare PowderThief - 1500mm – 20ml Tip

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 1500mm

Diameter of Outer Tube: 21mm

Nominal Tip Size: 20ml

Width of Handle: 103mm

Nominal Weight of Sampler: 238g

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 20

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21