

1. Part Number: 8063H-300S

2. Photograph:



3. Product Information:

Description:	SteriWare® Powder Lance - 300mm – Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	401mm
Blade Length:	300mm
Diameter of Blade:	22mm
Nominal Weight of Sampler:	63g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. Change Control: Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin: Made and packed in the UK.



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21