

1. Part Number: 8063H-500

2. Photograph:



3. Product Information:

Description:	SteriWare® Powder Lance - 500mm
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	601mm
Blade Length:	500mm
Diameter of Blade:	22mm
Nominal Weight of Sampler:	80g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. Change Control: Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21