

**1. Part Number:** 8063H-1000S

**2. Photograph:**



**3. Product Information:**

Description: SteriWare® Powder Lance - 1000mm - Sterile

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 953mm

Blade Length: 851mm

Diameter of Blade: 22mm

Nominal Weight of Sampler: 120g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

**4. Change Control:** Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:** Made and packed in the UK.

## 6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

## 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21