

**1. Part Number:** A316-500**2. Photo:****3. Product Information:**

Description	PharmaScoop® 500ml
Nominal Volume:	500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	240mm
Body Length:	140mm
Body Diameter:	95mm
Nominal Weight:	346g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20