

1. Part Number: A316-2500

2. Photo:



3. Product Information:

Description	PharmaScoop® 2500ml
Nominal Volume:	2500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	320mm
Body Length:	215mm
Body Diameter:	185mm
Nominal Weight:	964g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20