

**1. Part Number:** A316-5000

**2. Photo:**



### 3. Product Information:

Description	PharmaScoop® 5000ml
Nominal Volume:	5000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	355mm
Body Length:	225mm
Body Diameter:	220mm
Nominal Weight:	1165g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20