


1. Part Number: 8045R-2500S

2. Photograph:



3. Product Information:

| | |
|---|--|
| Description: | SteriWare PharmaScoop® RED 2500ml - sterile |
| Use: | Single use, disposable scoop |
| Material of Construction: | High Impact Polystyrene (HIPS) – rigid, red <ul style="list-style-type: none"> • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 2500ml |
| Overall Length: | 357mm |
| Handle Length: | 147mm |
| Maximum Width: | 166mm |
| Nominal Weight: | 224g |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Number of Scoops per Box: | 60 (4 boxes of 15 scoops) |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX rated |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-19