

Cole-Parmer®

Cole-Parmer® Flask Shakers - WA-200

- Vigorous shaking action
- Timed operation or continuous running
- Robust construction
- Ideal for extractions



Cole-Parmer® Flask Shakers - WA-200



Valuable time can be taken up mixing bottles and flasks by hand. Let our flask shaker take the strain. Holding up to eight flasks or bottles, up to 500ml capacity, it creates a vigorous mixing action by simulating hand shaking. It is especially useful for applications where prolonged shaking is required as it can be left on continuously and won't get tired!

Robust construction, mounted on four rubber feet to absorb vibration and prevent unnecessary movement on the bench. Analogue timer covers 10 to 60 minutes with a manual override. Electronic feed-back control ensures a constant speed irrespective of load.

Supplied with two side arms, eight clamps and Allen key. An optional extension kit is available which consists of two extension arms and four extra clamps increasing the capacity to twelve bottles or flasks.

Technical Specification

Specification	WA-200
Speed range	8 to 800 oscillations/min
Amplitude, mm	1.5
Dimensions, mm (w x d x h)	780 x 270 x 240
Maximum load, kg	3
Net weight, kg	8.3
Electrical supply	230V, 50Hz, 50W
IP Rating	31

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Flask shaker with two side-arms, eight clamps	51900-02	WA-200	WA-200	SF1
Spare clamp	99955-87	WA-200	N/A	SF1/1
2 x accessory side - arms and 4 x clamps. Increases the capacity to twelve bottles or flasks	99955-88	WA-200	N/A	SF1/2

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials
an antylia scientific company

Find out more at antylia.com