

Wilford Industrial Estate **Ruddington Lane** Wilford Nottingham **NG11 7EP**

> Tel: 0115 982 1111 Fax: 0115 982 5275

TECHNICAL DATA SHEET SLE7572

Date of issue: 31 August 2019

Product Code: SLE7572

Description



150ml container, in polystyrene, patient labelled with a metal screw cap containing an integral flowed in seal. Product supplied in 120's, bagged in 40's.

Declaration of Conformity

This product comes under the remit of the In Vitro Medical Devices Directive 98/79/EC and is classed as a general IVD product under annex III (specimen receptacle as per article 1(b) of the directive).

As the manufacturer we hereby declare that the requirements and obligations, as described in annex III sections 2 to 5 of the Directive have been met.

Authorised by: Andrew Littlewood

QADOCPROD: SLE7572 PAGE 1 of 2







Microbiological Status	Manufactured using an aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
Material of construction (body)	Polystyrene
Cap Liner Material	Daraseal
Material of construction (Cap shell)	Tinplate
Volume (ml)	150 ml
Bottle Dimensions (mm)	114.8 x 49 mm (approx.)
Operational temperature Range	-10 to 60 Degrees Centigrade
Inner Pack Quantity	40 in re sealable gripper bag
Outer Box Quantity	120

QADOCPROD: SLE7572 PAGE 2 of 2, ISSUE DATE: 31 August 2019

Authorised by: Andrew Littlewood



