

Wilford Industrial Estate **Ruddington Lane** Wilford Nottingham **NG11 7EP**

> Tel: 0115 982 1111 Fax: 0115 982 5275

18 July 2019

TECHNICAL DATA SHEET SLE7584

Date of issue: 18 July 2019

Product Code: SLF7584

Description



250ml container, with Printed Label, in polystyrene, with a metal screw cap with integral flowed in seal. Product supplied in 48's, bagged in 24's. CE Marked. Class I IVD.

Declaration of Conformity

This product comes under the remit of the In Vitro Medical Devices Directive 98/79/EC and is classed as a general IVD product under annex III (specimen receptacle as per article 1(b) of the directive).


As the manufacturer we hereby declare that the requirements and obligations, as described in annex III sections 2 to 5 of the Directive have been met.

QADOCPROD: SLE7584 PAGE 1 of 2 Authorised by: Andrew Littlewood









Microbiological Status	Manufactured using an aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
Material of construction (body)	Polystyrene
Cap Liner Material	Daraseal
Material of construction (Cap shell)	Tinplate
Volume (ml)	250ml
Bottle Dimensions (mm)	122 x 60mm (approx)
Operational temperature Range	-10 to 60 Degrees Centigrade
Inner Pack Quantity	24 in re sealable gripper bag
Outer Box Quantity	48

QADOCPROD: SLE7584 PAGE 2 of 2, ISSUE DATE 18.07.2019

Authorised by: Andrew Littlewood



