

Material description: BUFFER SOLUTION pH 2  
CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE1372	BUFFER SOLUTION pH 2 1L.
	CHE1374	BUFFER SOLUTION pH 2 5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
pH @ 10°C	pH			1.98	
pH @ 15°C	pH			2.02	
pH @ 20°C	pH	1.99	2.01	2.00	
pH @ 25°C	pH			2.00	
pH @ 30°C	pH			2.00	
pH @ 35°C	pH			2.02	
pH @ 40°C	pH			2.01	
pH @ 45°C	pH			2.02	
pH @ 50°C	pH			2.02	
Traceable to NIST					

Revision date: 26/11/2021

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2634/MS/20211126