

Material description: BUFFER SOLUTION pH 4

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE1376	BUFFER SOLUTION pH 4 1L.
	CHE1378	BUFFER SOLUTION pH 4 5L.
	CHE3962	BUFFER SOLUTION pH 4 500ml.
	CHE4126	BUFFER SOLUTION pH 4 20L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
pH @ 10°C	pH			4.00	
pH @ 15°C	pH			4.00	
pH @ 20°C	pH	3.99	4.01	4.00	
pH @ 25°C	pH			4.01	
pH @ 30°C	pH			4.01	
pH @ 35°C	pH			4.02	
pH @ 40°C	pH			4.03	
pH @ 45°C	pH			4.04	
pH @ 50°C	pH			4.05	
Traceable to NIST					

Revision date: 30/09/2021

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0119/MS/20210930