



Material description: EOSIN Y 24050 (C.I. 45380-water & alcohol soluble)

CAS No: 17372-87-1

Supplied as: Product code Description
CHE1916 EOSIN Y 24050 (C.I. 45380-water & alcohol soluble) 25g.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 80 | 88 | | |
| Wavelength of max absorption | nm | 515 | 517 | | |
| Specific Absorption | | 1120 | 1235 | | |
| Loss on Drying | % | | 8 | | |

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3105/MS/19890101