



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: GLYCINE pure  
CAS No: 56-40-6

Supplied as:	Product code	Description
	CHE2082	GLYCINE pure 100g.
	CHE2084	GLYCINE pure 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98.5			
Chloride	ppm		70		
Sulphate	ppm		65		
Heavy Metals; as Pb	ppm		20		
Loss on drying	%		0.2		

Revision date: 26/09/2024

Issue date: 18/09/2025

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0037/MS/20240926