



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: HYDROCHLORIC ACID 32% w/w tech.

CAS No: 7647-01-0

Supplied as: Product code Description
CHE2132 HYDROCHLORIC ACID 32% w/w tech. 5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	31	33		
Wt. per ml.	g	1.155	1.165		

Revision date: 28/06/2005

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3106/MS/20050628