



Material description: PHENOLPHTHALEIN SOLUTION (for milk testing)

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2852	PHENOLPHTHALEIN SOLUTION (for milk testing) 500ml.
	CHE2854	PHENOLPHTHALEIN SOLUTION (for milk testing) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Clear colourless to Pink					pH 8.3-10.0

Revision date: 12/05/2020

Issue date: 18/09/2025

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1627/MS/20200512