

Material description: POLYETHYLENEGLYCOL 400

CAS No: 25322-68-3

Supplied as: Product code
CHE2884

Description
POLYETHYLENEGLYCOL 400 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Average molecular weight		380	440		
Water content	% w/w		1.0		
Colour	Hazen		50		
pH 5% aqueous solution	pH	4.5	7.5		

Revision date: 11/06/2020

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1880/MS/20200611