

Material description: POTASSIUM BROMIDE pure
CAS No: 7758-02-3

Supplied as: Product code Description
CHE2896 POTASSIUM BROMIDE pure 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Bromate	%		0.001		
Chloride	%		0.25		
Iodide	%		0.01		
Moisture	%		0.5		

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0529/MS/19890101