

Material description: POTASSIUM CHLORIDE pure  
CAS No: 7447-40-7

<i>Supplied as:</i>	Product code	Description
	CHE2916	POTASSIUM CHLORIDE pure 3Kg.
	CHE2918	POTASSIUM CHLORIDE pure 500g.
	CHE2920	POTASSIUM CHLORIDE pure 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Lead	%		0.001		
Sulphate	%		0.03		
Moisture	%		1		

Revision date: 03/11/1997

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0531/MS/19971103