



Material description: SULPHURIC ACID 90% w/w (for milk testing)
CAS No: 7664-93-9

Supplied as: Product code Description
CHE3708 SULPHURIC ACID 90% w/w (for milk testing) 2.5L.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-------------|-------|---------------|---------------|---------------|-------|
| Strength | % w/w | 89.5 | 91.5 | | |
| Wt. per ml. | g | 1.812 | 1.818 | | |
| Chloride | % | | 0.00002 | | |
| Nitrate | % | | 0.00002 | | |

Revision date: 25/03/2009

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2400/MS/20090325