



Material description: FORMALDEHYDE SOLUTION 40% w/v pure
CAS No: 50-00-0

Supplied as: Product code Description
CHE2036 FORMALDEHYDE SOLUTION 40% w/v pure 2.5L.
CHE5032 FORMALDEHYDE SOLUTION 40% w/v pure 25L.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-----------------|-------|---------------|---------------|---------------|-------|
| Strength | % w/v | 39 | 41 | | |
| Wt. per ml. | g | | | 1.08 | |
| Iron | % | | 0.002 | | |
| Acidity; formic | % | | 0.05 | | |

Revision date: 28/06/1990

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0332/MS/19900628