

# Pharmacy Refrigerators

— **HYC-68/68A**



## Ergonomic Design



### Audible and Visual Alarms

Audible buzzer and visible flashing light



### Accurate Temperature Control



### Multiple Alarm Functions

High/low temperature, power failure, door ajar, sensor error, low battery, and remote alarm



### Interior LED



### Air Cooling, Auto Defrost

Superior cabinet temperature uniformity, auto defrost



### Lockable Design

## Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
Tel: +86-0532-88935593  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint



Haier Biomedical  
International



Haier Biomedical  
International

## Specifications



### Technical Data

Cabinet Type	Upright
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise (dB (A))	41

### Performance

Temp Range (°C)	2~8
-----------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply (V/Hz)	220~240/50	115/60
Power (W)	140	150
Electrical Current (A)	0.85	1.7

### Dimensions

Capacity (L/Cu.Ft)	68/2.4
Net/Gross Weight (approx)	38/40 (kg) 83.8/88.2 (lbs)
Interior Dimensions (W*D*H)	415*385*505 (mm) 16.3*15.2*19.9 (in)
Exterior Dimensions (W*D*H)	495*580*660 (mm) 19.5*22.8*26.0 (in)
Packing Dimensions (W*D*H)	540*640*710 (mm) 21.3*25.2*28.0 (in)
Container Load (20'/40'/40'H)	102/204/204



HYC-68



HYC-68A

### Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote Alarm	Y

### Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves /Drawers	2/1

### Other

Certification	CE	UL
---------------	----	----