

SLS Lab Pro Stirred Water Baths

SLS Lab Pro Stirred Water Baths

The SLS Lab Pro Range of Stirred Water Baths 8, 14, 22 and 56L provide class leading performance and reliability for scientific, laboratory and industrial applications for consistent temperature control.

- Precision PID temperature control, increments 0.1°C
- Bold LED temperature display
- Stainless steel tank
- Featuring a rugged compact sized circulator; to the rear a blanking hole accommodates immersion coolers
- Circulator can easily be removed for cleaning by undoing two thumb wheel securing screws
- Painted surfaces feature antibacterial finish
- Safety features, low liquid detection, with FILL messaging
- Factory temperature calibration, fully tested
- Accessories available, enquire for details



Environment, manufacturing facility

The range is designed and manufactured in the UK minimising 'transportation and distribution' miles. Wherever possible materials re-cycled or re-used. We aim to minimise energy use and investment in low energy plant, processes and equipment.

Standard specification	
Construction	Stainless steel 304 spec tank. Painted surfaces feature "anti-bacterial technology".
Temperature range, °C	Ambient +5 to 99, PID digital LED display showing actual or set temperature, resolution 0.1°C. When used in conjunction with accessory immersion cooler provides cooling from ambient to -15°C
Stability, °C	±0.5
Uniformity, °C	±0.01
Safety	Low liquid, "FiLL" message in display
Voltage, V	230
Warranty, years	3

	Overall Dims w x d x h, mm	Tank Dims w x d x h, mm	Tank max capacity, L	Energy use, kWh	Packed size, cm/Weight, kg w x d x h, cm/kg
SLS1170	271 x 392 x 170	129 x 300 x 150	8	0.5 @ 37°C	45 x 41 x 32/9.2
SLS1172	361 x 332 x 170	219 x 298 x 150	14	0.5 @ 37°C	45 x 41 x 32/10.7
SLS1174	537 x 332 x 170	395 x 298 x 150	22	0.5 @ 37°C	62 x 41 x 37/12.7
SLS1177	660 x 540 x 290	490 x 500 x 200	56	1.3 @ 37°C	69 x 47 x 79/23

Lift Off Stainless Steel Lids (not included)	
Code	Description
SLS1186	Lid for use with 8L baths
SLS1188	Lid for use with 14L baths
SLS1190	Lid for use with 22L baths
SLS0998	Lid for use with 56L baths

00229