

SLS Lab Basics Pipettes				SLS Lab Pro Pipette Tips					
SKU	Channel	Volume Range, μ L	Compatible Tips	Racked			Reload	Bulk	
				Sterile, Filter	Sterile	Non Sterile			
SLS4900	Single	0.1 - 2.5	10 μ L, 10XL	SLS5146, SLS5148	SLS5112, SLS5110	SLS5124, SLS5122	SLS5136, SLS5134	SLS5234, SLS5232	
SLS4902	Single	0.5 - 10	10 μ L, 10XL	SLS5146, SLS5148	SLS5112, SLS5110	SLS5124, SLS5122	SLS5136, SLS5134	SLS5234, SLS5232	
SLS4910	Single	2 - 20	10XL	SLS5148	SLS5110	SLS5122	SLS5134	SLS5232	
SLS4904	Single	5 - 50	20 μ L, 200 μ L*	SLS5150	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4906	Single	10 - 100	100 μ L, 200 μ L*	SLS5152	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4912	Single	20 - 200	200 μ L, 300 μ L	SLS5154, SLS5156	SLS5114, SLS5116	SLS5126, SLS5128	SLS5138, SLS5140	SLS5236, SLS5238	
SLS4908	Single	100 - 1000	1000 μ L, 1250 μ L	SLS5158, SLS5160	SLS5118, SLS5120	SLS5130, SLS5132	SLS5142, SLS5144	SLS5240, SLS5242	
SLS4914	Single	500 - 5000	-	-	-	-	-	SLS4042**	
SLS4916	Single	1000 - 10000	-	-	-	-	-	SLS4044**	
SLS4926	8	0.5 - 10	10 μ L, 10XL	SLS5146, SLS5148	SLS5112, SLS5110	SLS5124, SLS5122	SLS5136, SLS5134	SLS5234, SLS5232	
SLS4918	8	5 - 50	20 μ L, 200 μ L*	SLS5150	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4920	8	10 - 100	100 μ L, 200 μ L*	SLS5152	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4922	8	20 - 200	200 μ L, 300 μ L	SLS5154, SLS5156	SLS5114, SLS5116	SLS5126, SLS5128	SLS5138, SLS5140	SLS5236, SLS5238	
SLS4924	8	30 - 300	200 μ L, 300 μ L	SLS5154, SLS5156	SLS5114, SLS5116	SLS5126, SLS5128	SLS5138, SLS5140	SLS5236, SLS5238	
SLS4936	12	0.5 - 10	10 μ L, 10XL	SLS5146, SLS5148	SLS5112, SLS5110	SLS5124, SLS5122	SLS5136, SLS5134	SLS5234, SLS5232	
SLS4928	12	5 - 50	20 μ L, 200 μ L*	SLS5150	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4930	12	10 - 100	100 μ L, 200 μ L*	SLS5152	SLS5114	SLS5126	SLS5138	SLS5236	
SLS4932	12	20 - 200	200 μ L, 300 μ L	SLS5154, SLS5156	SLS5114, SLS5116	SLS5126, SLS5128	SLS5138, SLS5140	SLS5236, SLS5238	
SLS4934	12	30 - 300	200 μ L, 300 μ L	SLS5154, SLS5156	SLS5114, SLS5116	SLS5126, SLS5128	SLS5138, SLS5140	SLS5236, SLS5238	

* Non-filter tips only

* SLS Select bulk pipette tips