Passive Accessories

73 Series Accessories

Jenway also offer an extensive range of passive accessories to use with the 73 series spectrophotometers. For cuvettes with a variable path length between 10 and 100mm we offer an adjustable path length cuvette holder. **Part code: 630 005**

For small sample volumes down to $50\mu l$ we offer a micro-cuvette holder. Part code: $630\ 304$

For large sample volumes we offer a test tube holder which can hold test tubes up to 100mm in height with diameters of 13, 16 and 24mm Part code: 637 071. Jenway also offer a boiling tube holder which can hold boiling tubes up to 150mm in height with a diameter of 25mm. Part code: 735 601

Printers

For instant results Jenway offer a printer which can be easily fitted into the top of the spectrophotometer to save bench space. The integrated printer enables spectrum and kinetics curves to be printed as well as results tables and analysis summary. The user name and method parameters are printed on all results.

Part code: 660 101

Ordering Information

Part Code	Description
660 101	Internal printer (packed)
735 401	Automatic 8 position cuvette holder
735 201	Sipper pump
735 301	Peltier
735 701	Combined sipper/Peltier
630 204	10x10mm path length cuvette holder
630 005	10 to 100mm path length cuvette holder
637 071	16/24mm test tube holder
630 304	Micro-cuvette holder with reduced aperture
736 201	Water heated 10x10mm path length cuvette holder
735 601	Boiling tube holder
735 001	Dust cover
019 146	4GB USB Memory stick for external memory
	(7310/7315 only)
012 075	Tungsten halogen lamp (7300/7310)
012 146	Xenon lamp module (7305/7315)







