

250 μL, Model 725 RN SYR, Large Removable NDL, 22s ga, 2 in, point style 2

Specification Sheet (Part/REF # 80730)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

Product Specifications

Volume	250 μL
Syringe Type	Microliter
Syringe Series	700
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F 50°C / 122°F
Minimum Temperature	10°C/50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)



Visit $\underline{www.hamiltoncompany.com/support}$ to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055



Specifications are subject to change without notice	Spec. Version B
Sales Quantity	1 syringe per pack
Std Needle Dead Volume	1.13 μL
Needle Outer Diameter	0.0283 inches (0.72 mm)
Needle Inner Diameter	0.0059 inches (0.15 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)

https://www.hamiltoncompany.com/laboratory-products/syringes/80730



US: +1-775-525-3000 EU: +40-356-635-055