

50 μL, Model 1705 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2

Specification Sheet (Part/REF # 80930)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

Product Specifications

Volume	50 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C/50°F
Barrel Inner Diameter	0.041 inches (1.03 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)



Visit $\underline{www.hamiltoncompany.com/support}$ to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055



Needle Outer Diameter 0.0283 inches (0.72 mm) Std Needle Dead Volume 1.13 μL Sales Quantity 1 syringe per pack	ersion B
Needle Outer Diameter 0.0283 inches (0.72 mm)	
0.0000 11101103 (0.10 11111)	
Needle Inner Diameter 0.0059 inches (0.15 mm)	

HAMILT®N°

US: +1-775-525-3000 EU: +40-356-635-055