

## 25 μL, Model 1702 LT SYR, NDL Sold Separately

Specification Sheet (Part/REF # 80201)



The 1700 series syringes are the Gastight version of the original 700 series syringes. This series of syringes is ideal for dispensing volumes from 1  $\mu$ L up to 500  $\mu$ L. A variety of different terminations and needle options are available in this series.

## **Product Specifications**

Volume	25 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	Luer Tip (LT)
Needle	Kel-F Hub
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes to 115 °C
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C/50°F
Barrel Inner Diameter	0.029 inches (0.729 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)



Visit  $\underline{www.hamiltoncompany.com/support}$  to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055



Std Needle Dead Volume	N/A	
Sales Quantity	1 syringe per pack	
Specifications are subject to change without notice		Spec. Version B

https://www.hamiltoncompany.com/laboratory-products/syringes/80201



Visit <u>www.hamiltoncompany.com/support</u> to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055