

# 26s gauge, Large Hub RN Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 7779-02)



Large Hub Removable (RN) needles can be used with removable (RN)-style syringes that are 250  $\mu$ L and larger. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

|                       |                    |
|-----------------------|--------------------|
| <b>Gauge</b>          | 26s gauge          |
| <b>Point Style</b>    | 2                  |
| <b>Needle Length</b>  | 2 inch (51 mm)     |
| <b>Autoclavable</b>   | Yes                |
| <b>Sales Quantity</b> | 6 needles per pack |
| <b>Hub</b>            | Large RN           |

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/needles/7779-02>